KNGMED

KNGMED MEDİKAL LTD. ŞTİ.
29 Ekim Mah. 10007 Sk.
No:26/B, Menemen, İZMİR,
TÜRKİYE
Tel:0(232)8334262

www.kngmed.com
info@kngmed.com&kng.info@yahoo.com

KNGSORB*Sodalime Carbon Dioxide Absorbent*

LOT

REF

UDI

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8 699947 601207

Intended Use:

KNGSORB in the close circuits of anesthesia systems, used to absorb CO₂.

Principle of Operation:

KNGSORB binds the CO₂ contained in the breathing gas in an exothermic reaction. The hydroxides are converted into carbonates and water.

Instruction for Use:

The canister must be filled with KNGSORB appropriately, when placed to the machine, it will be ready to use.

Replacing:

KNGSORB Sodalime should be replaced when at least one of the following occurs:

- When the inspiratory CO₂ concentration reaches approximately 5mmHg,

- When about 2/3 of the soda lime turns purple,

- When sodalime stays in the anesthesia device for 1 week, it must be changed even if it has not been used.

Warnings and Precautions:

The product should be kept in its original packaging.

-The packaging lid should always be kept closed to prevent exposure to moisture.

-Irritating in contact with skin.

-Suitable protective gloves should be used to prevent skin exposure.

-This product must be kept away from children.

-The product is for single use only. It should be replaced with a new one as soon as it turns purple.

-Check the sodalime status and machine indicators frequently.

-Products with damaged packaging should not be used.

-Dispose of as medical waste in accordance with local legislation.

-Inhaled anesthetic agents in contact with strong bases have been reported to degrade to produce breakdown products of unknown toxicity in humans. For further information, contact the manufacturer of the inhaled anesthetic agent in question.

-Technicians should be warned to seek medical advice in case of skin irritation.

-The product must not be used by lay persons.

-Do not use expired products.

-The product must not be stored under inadequate storage conditions.

-The canister must be sufficiently filled with Sodalime.

-The Sodalime medical device used in anesthesia devices must be used in a clean and hygienically safe environment.

EN