

reusable silicon anesthesia mask



*reusable silicon
anesthesia mask*



Raw material

- › 100% anti-allergenic, odorless, medical silicone

Design

- › Leakage preventive inward-oriented grooved design
- › Transparent form

Features

- › Universal 22 mm tube

Usage

- › Reusable
- › All models can be autoclaved repeatedly at 134° C
- › Suitable for all patient circuits

Size

- › Five different sizes (No. 1 thru 5)

1. No:00 item no: #7600125-00
2. No:01 item no: #7600125-01
3. No:02 item no: #7600125-02
4. No:03 item no: #7600125-03
5. No:04 item no: #7600125-04
6. No:05 item no: #7600125-05



Tel : +90 232 833 42 62

Fax: +90 232 833 42 63

info@kngmed.com

kng.info@yahoo.com

10001 Sokak No.28 / A - 10

Çınar Sanayi Sitesi

Ulukent / Menemen / İzmir / TÜRKİYE