

broncho alveolar mucus extractor

Mucus Extractor is utilized to obtain sputum and/or tracheal aspirate specimens to identify possible bacteria in the lungs. Allows for easy, convenient collection of mucous specimen samples. After specimen collection, the two connectors are joined to close the system.



40cc



15cc

Design

- › Special design

Features

- › Male/female connection port
- › Sterile single package
- › Scale shows sample amount
- › Name tag

Usage

- › Single use only
- › Suitable for all bronchoscope devices

Size

- › 40cc

Design

- › Special design

Features

- › Ch/Fr 14 and Ch/Fr 16 connectors
- › Sterile single package
- › Scale shows sample amount
- › Name tag

Usage

- › Single use only
- › Suitable for all bronchoscope devices

Size

- › 15cc

15cc item no: #7600057-15
40cc item no: #7600057-40



Tel : +90 232 833 42 62
Fax: +90 232 833 42 63

info@kngmed.com
kng.info@yahoo.com

10001 Sokak No.28 / A - 10
Çınar Sanayi Sitesi
Ulukent / Menemen / İzmir / TÜRKİYE

broncho alveolar mucus extractor