

KNGMED





KNGMED MEDİKAL was established in Ankara in 2008 and initially joined medical sector as an importer company.

In order to start production of medical devices; moved to İzmir in 2009 and started production of sleep laboratory consumables at its first production site in Konak.


After moving to second production site of 600m² in 2010, CPAP Masks and Co₂ Absorbent were added to product range, and then Internal Nasal Splints were produced in same year.

Moved to its current location at the beginning of 2019 and still produces with 30 personnel in factory on an area of 3750m². Serves not only Medical Devices industry, but also the cosmetics industry.

In addition, our company, which has dealers in 81 cities in Turkey, exports to 35 countries abroad.

KNGMED Medical product range consists of Soda Lime, CPAP- Bi-PAP Masks, Ventilator Masks and Reusable Silicone Anesthesia Masks, Internal Nasal Splints, Silicone Menstrual Cups, CPAP – Bi-PAP Circuits.

Our company, which does not compromise on quality, safety and comfort of use in production, aims to be a close follower of rapidly technology today.



**Oro-Nasal
Masks**



**Total Face
Masks**




Nasal Masks



**Reusable Silicone
Anesthesia Masks**



Accessories



**CPAP BI-PAP
Circuits**



**Silicone Internal
Nasal Splints**




**Menstrual
Cups**



**Sodalime
Co2 Absorber**

INOVATION

new model



One MASK ADAPTABLE uses



Comfort you will feel

- Medical grade silicone sealing cushion provides optimal comfort for the patients.
- Superlight weight and clear view to ease anxiety and claustrophobia.
- Slide-in elips for headstrap design for easy-to-use.
- Anti-Asphyxia Elbow is for patient safety; optional Standard Elbow good for ventilator application.
- Available sizes XS / S / M / L / XL

TIP!!!

To Choose mask

AAV* elbow; has an anti- asphyxia valve (AAV) for dedicated non-invasiventilators (use with single-limb circuits).

AAV elbow with leak; has an anti- asphyxia valve for dedicated non-invasive ventilators and intentional leak ports for sub-acute and homecare systems when intentional leak is needed.

Standard elbow (SE); is typically used with a dual limb circuit on ventilors equipped with integrated safety valve systems.



Double
Layer



Code	Description
7601894-XS	AAV* Elbow , intentional leak (Exhalation), Headgear, XSmall
7601894-S	AAV* Elbow , intentional leak (Exhalation), Headgear, Small
7601894-M	AAV* Elbow , intentional leak (Exhalation), Headgear, Medium
7601894-L	AAV* Elbow , intentional leak (Exhalation), Headgear, Large
7601894-XL	AAV* Elbow , intentional leak (Exhalation), Headgear, Xlarge



Code	Description
7601573-XS	AAV* Elbow ,Headgear, XSmall
7601573-S	AAV* Elbow Headgear, Small
7601573-M	AAV* Elbow ,Headgear, Medium
7601573-L	AAV* Elbow Headgear, Large
7601573-XL	AAV* Elbow ,Headgear, Xlarge
7601574-XS	Standard Elbow, Headgear,XSmall
7601574-S	Standard Elbow, Headgear,Small
7601574-M	Standard Elbow, Headgear,Medium
7601574-L	Standard Elbow, Headgear,Large
7601574-XL	Standard Elbow, Headgear,XLarge




Can adapt one mask to multiple ventilators with **interchangeable** elbows.

Code	Description
7601575-XS	AAV* Elbow , intentional leak (Exhalation), Headgear, XSmall
7601575-S	AAV* Elbow , intentional leak (Exhalation), Headgear, Small
7601575-M	AAV* Elbow , intentional leak (Exhalation), Headgear, Medium
7601575-L	AAV* Elbow , intentional leak (Exhalation), Headgear, Large
7601575-XL	AAV* Elbow , intentional leak (Exhalation), Headgear, Xlarge

555 Series

Oro Nasal Masks

Forehead Support



Code	Description
7600224-XS	AAV* Elbow , intentional leak (Exhalation), Headgear, XSmall
7600224-S	AAV* Elbow , intentional leak (Exhalation), Headgear, Small
7600224-M	AAV* Elbow , intentional leak (Exhalation), Headgear, Medium
7600224-L	AAV* Elbow , intentional leak (Exhalation), Headgear, Large
7600224-XL	AAV* Elbow , intentional leak (Exhalation), Headgear, Xlarge



Code	Description
7601276-XS	AAV* Elbow ,Headgear, XSmall
7601276-S	AAV* Elbow Headgear, Small
7601276-M	AAV* Elbow ,Headgear, Medium
7601276-L	AAV* Elbow Headgear, Large
7601276-XL	AAV* Elbow ,Headgear, Xlarge
7601277-XS	Standard Elbow, Headgear,XSmall
7601277-S	Standard Elbow, Headgear,Small
7601277-M	Standard Elbow, Headgear,Medium
7601277-L	Standard Elbow, Headgear,Large
7601277-XL	Standard Elbow, Headgear,XLarge



Can adapt one mask to multiple ventilators with **interchangeable** elbows.

Code	Description
7601278-XS	AAV* Elbow , intentional leak (Exhalation), Headgear, XSmall
7601278-S	AAV* Elbow , intentional leak (Exhalation), Headgear, Small
7601278-M	AAV* Elbow , intentional leak (Exhalation), Headgear, Medium
7601278-L	AAV* Elbow , intentional leak (Exhalation), Headgear, Large
7601278-XL	AAV* Elbow , intentional leak (Exhalation), Headgear, Xlarge

Oro Nasal Masks

Chin&Forehead Support


444 Series



Code	Description
7600904-S	Standard Elbow, Headgear,Small
7600904-M	Standard Elbow, Headgear,Medium
7600904-L	Standard Elbow, Headgear,Large



Code	Description
7600903-S	AAV* Elbow Headgear, Small
7600903-M	AAV* Elbow ,Headgear, Medium
7600903-L	AAV* Elbow Headgear, Large



Exhalation



Code	Description
7600905-S	AAV* Elbow , intentional leak (Exhalation), Headgear, Small
7600905-M	AAV* Elbow , intentional leak (Exhalation), Headgear, Medium
7600905-L	AAV* Elbow , intentional leak (Exhalation), Headgear, Large

333 Series

Oro Nasal Masks
 Chin Support




Code	Description
7600751-S	Standard Elbow, Headgear, Small
7600751-M	Standard Elbow, Headgear, Medium
7600751-L	Standard Elbow, Headgear, Large



Code	Description
7600750-S	AAV* Elbow Headgear, Small
7600750-M	AAV* Elbow ,Headgear, Medium
7600750-L	AAV* Elbow Headgear, Large



Exhalation



Code	Description
7600752-S	AAV* Elbow , intentional leak (Exhalation), Headgear, Small
7600752-M	AAV* Elbow , intentional leak (Exhalation), Headgear, Medium
7600752-L	AAV* Elbow , intentional leak (Exhalation), Headgear, Large

Oro Nasal Masks

Forehead Supported


222 Series



Code	Description
7600041-S	Standard Elbow, Headgear, Small
7600041-M	Standard Elbow, Headgear, Medium
7600041-L	Standard Elbow, Headgear, Large



Code	Description
7600040-S	AAV* Elbow Headgear, Small
7600040-M	AAV* Elbow ,Headgear, Medium
7600040-L	AAV* Elbow Headgear, Large



Exhalation

Code	Description
7600042-S	AAV* Elbow , intentional leak (Exhalation), Headgear, Small
7600042-M	AAV* Elbow , intentional leak (Exhalation), Headgear, Medium
7600042-L	AAV* Elbow , intentional leak (Exhalation), Headgear, Large

111 Series


Oro Nasal Masks



Code	Description
7600607-S	Standard Elbow, Headgear, Small
7600607-M	Standard Elbow, Headgear, Medium
7600607-L	Standard Elbow, Headgear, Large



Code	Description
7600606-S	AAV* Elbow Headgear, Small
7600606-M	AAV* Elbow ,Headgear, Medium
7600606-L	AAV* Elbow Headgear, Large



Exhalation

Code	Description
7600608-S	AAV* Elbow , intentional leak (Exhalation), Headgear, Small
7600608-M	AAV* Elbow , intentional leak (Exhalation), Headgear, Medium
7600608-L	AAV* Elbow , intentional leak (Exhalation), Headgear, Large

Accessories



FF-XS
Oro Nasal
Cushion
XSmall



FF-S
Oro Nasal
Cushion
Small




FF-M
Oro Nasal
Cushion
Medium



FF-L
Oro Nasal
Cushion
Large



FF-XL
Oro Nasal
Cushion
XLarge



KB-77 Headgear For Total Face Mask

KB-88 Headgear for Ora-Nasal and Nasal Masks


Accessories




TL-01
Standard Elbow



TL-02
AAV Elbow with Intentional Leak




TL-03
Standard Elbow with Bronchoscopy Port



TL-04
AAV Elbow

TL-05
Standard Elbow with Nebuliser Port



TL-06
AAV Elbow with Nebuliser Port



TC-01
Swivel replacement kit (22mm male-male)



EX-02
Disposable exhalation port with filter exhalation port with swivel (no filter)




EX-01
Disposable exhalation port with Oxygen enrichment attachment



EX-03
Disposable exhalation port with Oxygen enrichment attachment

Circuits



for Professional Use

CPAP BI-PAP Circuits

Code	Description
7601245	EVA

- Excellent clarity,
- Smooth bore for maximum airflow and reduced sound level,
- Odorless.



for Homecare Use

CPAP BI-PAP Hose

Code	Description
BD7601252	PE With TPE Cuff

- Strong, lightweight, and versatile hose for a variety of applications,
- Excellent clarity,
- Odorless,
- Smooth bore for maximum airflow and reduced sound level, a clear choice for disposable medical applications.

Breath



Code	Description
7601056-S	AAV* Elbow , intentional leak (Exhalation), Headgear, Small
7601056-M/L	AAV* Elbow , intentional leak (Exhalation), Headgear, Medium
7601056-XL	AAV* Elbow , intentional leak (Exhalation), Headgear, Large
7601055-S	Standard Elbow, Headgear,Small
7601055-M/L	Standard Elbow, Headgear,Medium/Large
7601055-XL	Standard Elbow, Headgear,XLarge
7601054-S	AAV* Elbow Headgear, Small
7601054-M/L	AAV* Elbow ,Headgear, Medium/Large
7601054-XL	AAV* Elbow Headgear, XLarge

Also can adapt interchangeable elbows. (See Accessories)

Can adapt one mask to multiple ventilators with **interchangeable** elbows.

Nasal Mask

Nap Series



Code	Description
7600217	With headgear, standard size

- SuperLight weight
- 360 rotating tubing connector
- Quick-release elbow
- 100% Medical Grade Silicone
- Flexible body
- Soft cushion
- Reconfigured Exhalation Vent Reduces Noise
- Secure Seal and Universal Sizing

Silicone Internal Nasal Splints



Code	Description
7600361-A	Sterile Silicone Internal Nasal Splint - Integral Airway

- Soft Medical Grade Silicone
- Funnel Shaped Airway design
- Designed to increase airflow
- Designed to minimize blockage
- Designed to support nasal septal structure
- Sterile-Ready-To-Use
- At health centers for professional use

Reusable Silicone Anesthesia Masks



Code	Description
760125-00	No:0
760125-01	No:1
760125-02	No:2
760125-03	No:3
760125-04	No:4
760125-05	No:5

- 100% Medical Grade Silicone
- Leakage preventive inward-oriented grooved design
- Transparent form
- Standard 22mm for sizes: 2,3,4,5
- Standard 15 for sizes: 0,1
- Reusable-Autoclavable
- Anatomical design

Sodalime

CO₂ Absorber



KNGSORB is a medical soda lime manufactured from mixture of Calcium and Sodium Hydroxides. It is supplied in the form of hard, porous regular rounded pellets which have been specifically processed to absorption CO₂ capacity. One kilogram of KNGSORB will absorb a minimum of 130 litres of CO₂.


Code	Description
7600064-5	5 Kg Barrels , White to Violet

With the irregular shape, minimises the amount of dust produced during transit allowing the product to reach you in perfect condition for use.

KNGSORB does not contain potassium Hydroxide or Barium Hydroxide and only low level of Sodium Hydroxide.

High CO₂ absorption capacity pellet from prevents channel formation, uniform utilization reducing your costs through longer life time little dust generation creating a safer environment for your personnel and protects reduction of undesirable reaction products when accidentally dried out.

Menstrual Cup





Ivno The menstrual cup is a new invention that will put an end to the use of pads and tampons. Ivno is easily inserted at the entrance of the vagina and collects the menstrual blood from the uterus in a way that does not leak out. It prevents all kinds of leakage with its special design in the shape of a small funnel. It is a hygienic, body and budget friendly product that can be used safely for many years.



Code	Description
7601269-S	Small
7601269-M	Medium
7601269-L	Large

KNGMED



www.kngmed.com - info@kngmed.com +90 232 8334262

29 Ekim Mh. 10007 Sk. No:26/B Menemen / İZMİR - TÜRKİYE